


Equipment ID: B4/CM/190



FirstSource Laboratory Solutions LLP
 Plot No.2, Block-D, Street No.3
 Kakatiyanagar, Hubsiguda, Hyderabad-500 007.
 Ph: 040-27175520-22 Fax: 040-27175525
 Email: sales@firstsourcels.com, Web: www.firstsourcels.com

INSTALLATION REPORT

CRM PM No :		Spinlab ID:	Engineer Name: Jordan	Date: 20.09.23
Account :	112	Customer ID:	Letter	Date:
Address :	Gachibowli	Principal: Eutech	Invoice no:	Date:
State :	TS	Product Name: 2700/CM	Job Completed: Yes <input checked="" type="checkbox"/> No <input type="checkbox"/>	
Pin Code:	500032	System User Details		
SID In-Charge Name:	E. VENKAT REDDY	Name	Designation	Mobile
Designation:	Sr. Manager	NA		
Department:	HH - Production			
Phone:	-			
Fax:	-			
mobile:	9705052221			
Email:	VenkatReddy@indirmune.com			
Website:	-			

Details : Installation 8/NI- 3129560
 Printer- Antech :- MAX36KA 28040004

Installation Procedure

- (i) Installed the CON 2800
- (ii) Performed 1,00 - all okay, handed over to users.
- (iii) Calibrated - in limit.
- (iv) Gave demo to the users

AC Line Voltage Reading : P-N = 230 V P-E = 225 V E-N = 1 V UPS Stabiliser

Further Action Required : Make permanent earth connection before start using the instrument.

This is to certify that the above mentioned work has been carried out and the system is working to our satisfaction.

Customer Comments : Nil

			Engineer Signature :
Date	Time IN	Time OUT	
20.09.23	11:55	16:30	
			Customer Signature : (With Seal)
			Date: 22-09-23
K. Vijender Reddy (K. VIJENDER REDDY)			