



TRIDENT EQUIPMENTS PVT. LTD.

A - 105,, Kailas Industrial Complex, Hiranandani Godrej Link Road, Vikhroli (W),
Mumbai - 79. Tel.: 2518 1704 / 05 / 06 Telefax : 2518 5463 E-mail : info@tridentequipments.com

Service Report

No. **63401**

Name of Customer : <u>Indicus Immunologicals</u> <u>limited.</u>	Product : <u>Thermo Scientific</u>
Address : <u>Rakshapuram.</u>	Model(s) : <u>Oxon Sta A211</u>
Contact Person/Dept : <u>Mr. Sushakar Reddy.</u>	Serial No.(s) : <u>23128</u>
PO No. : <u>451100455</u>	Installation <input type="checkbox"/> Warranty <input type="checkbox"/> Chargeable <input type="checkbox"/> Courtesy <input type="checkbox"/>
Telephone : <u>9985684000</u>	Repaired at - Customer Site <input checked="" type="checkbox"/> Service Centre <input type="checkbox"/>
Email : <u>P-Sushakarreddy@indicus.com</u>	AMC : Preventive Maintenance Visit No. <u>-</u>
	Breakdown Visit No. <u>-</u>
	Date : <u>22-09-23</u>

Nature of Problem :

* Installation of Oxon Sta A211, (BU/PH/1A3),
PH meter.

Action Taken :

- * Installed all Accessories as per Instruction manual.
- * Configured all Parameters as per Installations Qualifications.
- * Performed Calibration and Verification with PH Simulator.
- * Performed Calibration and verification with PH buffers
- * Observed all values are within limits.

Engineer's Comment :

* Installation of Oxon Sta A211 is Successfully completed.

Customer's Comment :

Installation completed & found satisfactory.

M.N. Unod Kumar 22-09-23
Service Engineer's Name & Signature :

P. Sushakar Reddy
Sushakar 22-09-23
Customer's Name, Signature & Stamp :

Encl. : 8976898568