



TRIDENT EQUIPMENTS PVT. LTD.

A - 105,, Kailas Industrial Complex, Hiranandani Godrej Link Road, Vikhroli (W),
Mumbai - 79. Tel.: 2518 1704 / 05 / 06 Telefax : 2518 5463 E-mail : info@tridentequipments.com

Service Report

No. 63412

Name of Customer : <u>Indian Technological Limited</u>	Product : <u>Thermo Scientific</u>
Address : <u>Pokshapuram.</u>	Model(s) : <u>Oxon star A211</u>
Contact Person/Dept : <u>Mr. Venkat</u>	Serial No.(s) : <u>X82444</u>
PO No. : <u>4511004658</u>	Installation <input checked="" type="checkbox"/> Warranty <input type="checkbox"/> Chargeable <input type="checkbox"/> Courtesy <input type="checkbox"/>
Telephone : <u>9705052221</u>	Repaired at - Customer Site <input checked="" type="checkbox"/> Service Centre <input type="checkbox"/>
Email : _____	AMC : Preventive Maintenance Visit No. <u>-</u>
	Breakdown Visit No. <u>-</u>
	Date : <u>03.10.23</u>

Nature of Problem :

* Installation of Oxon star A211 (B8/PH/112)

Action Taken :

- * Installed all Accessories as per Instruction manual.
- * Configured all Parameters as per Installation qualifications.
- * Performed Calibration and verification with pH Simulator
- * Performed Calibration and verification with pH buffers
- * Observed all values are within limits
- * Demo and training given to users.

Engineer's Comment :

* Installation of Oxon star A211 is Successfully Completed.

Customer's Comment : Installation completed as per IQCOP document & found satisfactory.

Service Engineer's Name & Signature : M.V.Vinod Kumar 03.10.23

Customer's Name, Signature & Stamp : EVR Reddy 03.10.23

Encl. : 8976898568.