




SERVICE REPORT

 Tecpharma SYSTEMS INC. MIG-315, Balaji Nagar, Kukatpally, Hyderabad - 500 072. Tel : 40 23151666, 23151077 E-mail : krishna@tecpharma.in	Date : 14/11/2019	S/R : 1129
	Engineer's Name : D. KRISHNA KRISHNA	
	Mobile No. : 909 209 2726	
Customer's Name :	Type of Services	
Address : Indian Immunologics Ltd.	1. Installation <input checked="" type="checkbox"/>	6. Paid Service <input type="checkbox"/>
Raksha Padam	2. Breakdown <input type="checkbox"/>	7. Calibration <input type="checkbox"/>
Gachibowli Post	3. AMC <input type="checkbox"/>	8. IQ/OQ/PQ <input checked="" type="checkbox"/>
Hyderabad.	4. Warranty <input checked="" type="checkbox"/>	9. Application <input type="checkbox"/>
	5. Inspection <input type="checkbox"/>	10. Goodwill <input type="checkbox"/>
Tel. No. : Fax No. :	Job Completed : Yes <input checked="" type="checkbox"/> No <input type="checkbox"/>	
E-mail :		
Instrument	Model No.	Sr. No. / INST ID
H BIO ULTRA LOW TEMPERATURE FREEZER	DWHTLS285	190130772
Problem As Reported by the Customer : freezer installation		
Service Done : H Bio ULT freezer installation completed. freezer qualification done. During PQ 72 hours temperature mapping done without load and without freezer is maintain set temperature within acceptable criteria.		
Note : criteria.		
Spare Parts Used / Required	To Be filled by Customer	
- NA -	<input checked="" type="checkbox"/> Instrument Installed & Found Working Satisfactorily <input type="checkbox"/> Prevention Maintenance Carried Out Satisfactorily <input type="checkbox"/> Fault Rectified & Instrument Working Satisfactorily <input type="checkbox"/> Please await our written approval for parts replacement <input type="checkbox"/> Comments (If Any)	
Engineer's Signature : 	Name : D. Krishna	
14/11/2019	Date : 14.11.19	Customer Signature
* Only for parts Changed & all other charges will be extra as applicable		

Note : Parts replaced are chargeable except during warranty period. Consumables used for the instrument are not covered by warranty & hence are chargeable. Parts replaced due to mishandling during operation or otherwise will also be charged in every case. In case of any difficulty or suggestion, please contact Service Manager.