

Installation Report			
User Name :	Mr. P. Satheesh	Date :	09/04/2021
Institution/Company :	Analytical Lab, IIL - Gachibowli, Hyderabad.	Mobile :	9133644007
Email :	Samuel@indimmune.com.		
Customer PO No. and Date :		Model No :	Digital Hotplate DAIHAN HP-LP
Invoice No and Date :		SI. No :	040198 ⁴ 2024001
Installation Date :	09/04/2021	Bar Code :	- NA -
Inspection :	Received And checked the Unit "found OK."		
Unit performance :	The unit was performing Well at 190°C.		
Guidelines / SOP :	1. Use gloves while handling the unit at High temperatures. 2. Clean the Surface of the unit after use. 3. Use operating manual for any operating procedure. 4. Keep the unit in flat Surface Area Always.		
Start Time :	End Time :	Total Time Spent :	
Local Accessories if any :			
Part No	Description	Qty	Remarks
	16A, 3 pin top plug	01	Replaced.
Remarks : Explained operational procedure to user			
SRICO Engineer :	Unit was checked at 100°C and 190°C, Works good. Unit installed Successfully.	User :	
Name and Signature of Engineer :	K. Sai Krishna + B. Krishna <i>K. Sai Krishna</i>	Name and Signature of User :	[SAMUEL] Samuel
Mobile :	7093001023	Company Stamp	
Email :	Sai.krishna @srico-labworld.com		

Scientific Research Instruments Company Private Limited

Telephone No : +91 9900055879

E-Mail : service@srico-labworld.com

Website : www.srico-labworld.com

Branches: Bangalore :

BBSR :

Delhi :

Hyderabad :

Thane :

Customer Copy :

Head Office Copy :

Branch Copy :

Service Engineer Copy :

Internal office use :