

Installation Report			
User Name : <u>Mr. Satish Kumar Manjunatha</u>	Date : <u>06-03-2020</u>		
Institution/Company : <u>Indian Immunologicals Limited.</u>			
<u>Gachibowli Hyderabad</u>	Mobile : <u>9700642870</u>		
Email: <u>m.Satishkumar@indimmune.com</u>			
Customer PO No. and Date : <u>451103104 Dt. 12-02-2020</u>	Model No : <u>2mag mix-4mg</u>		
Invoice No and Date : <u>1920/577 Dt. 24-02-2020</u>	Sl. No : <u>08357003/0517</u>		
Installation Date : <u>unit is in 06-03-2020</u>	Bar Code : _____		
Inspection : <u>unit is in Good condition no damage found</u>			
Unit performance : <u>after unpacking unit - and unit is in working</u>			
Guidelines / SOP : <u>as its specifications</u>			
<u>magnetic stirrer should be kept on flat surface. unit performance is checked and operation of it shown to the users</u>			
Start Time : <u>12:30 pm</u>	End Time : <u>2:00 pm</u> Total Time Spent : <u>1h30 mins</u>		
Local Accessories if any :			
<b>Part No</b>	<b>Description</b>	<b>Qty</b>	<b>Remarks</b>
	<u>SA Plug TOP</u>	<u>01</u>	<u>Replaced</u>
Remarks : <u>Installation was found satisfactory</u>			
SRICO Engineer : <u>Ac Power supply should not fluctuations use UPS supply is better.</u>		User : <u>Satish Kumar Manjunatha</u> <u>Sr. Manager, Process Development</u> <u>R&amp;D, Gachibowli, Hyderabad.</u>	
Name and Signature of Engineer : <u>Pawan + Seal</u>		Name and Signature of User : <u>[Signature]</u> <u>06.03.2020.</u>	
Mobile : <u>8296681123</u>			
Email : <u>Pawan @srico-labworld.com</u>		Company Stamp	

Scientific Research Instruments Company Private Limited

Telephone No : +91 9900055879

E-Mail : [service@srico-labworld.com](mailto:service@srico-labworld.com)

Website : [www.srico-labworld.com](http://www.srico-labworld.com)

Branches: Bangalore :

BBSR :

Delhi :

Hyderabad :

Thane :

Customer Copy :

Head Office Copy :

Branch Copy :

Service Engineer Copy :

Internal office use :