

Installation Report

User Name : MR. PHANI KUMAR Date : 27/11/2019
 Institution/Company : INDIAN IMMUNOLOGICALS LTD, R&D Lab
Rakshapuram, Gachibowli, Hyderabad Mobile : _____
 Email: _____
 Customer PO No. and Date : _____ Model No : DAIHAN MS-MPS
 Invoice No and Date : _____ SI. No : 11000462194 Foot
 Installation Date : 27/11/2019 Bar Code : -NA-

Inspection : Checked for Transport Damages, nothing found
 Unit performance : All stirring points and Remaining functions working
 Guidelines / SOP : (1) Install Instrument in ventilated area. (2) Do not place near any magnetic or metal materials. (3) Use Rubber mat always to Avoid slipping of flasks and also to avoid scratches. (4) use safe power supply. (5) Avoid splitting of chemicals on Display.

Start Time : _____ End Time : _____ Total Time Spent : _____

Local Accessories if any :

Part No	Description	Qty	Remarks
	<u>3 PIN PLUG TOP 6A</u>	<u>01</u>	<u>Replaced</u>

Remarks : Installation is Completed against above PO.

SRICO Engineer : Use Instrument User : _____
Between +5°C to +40°C
Only.
 Hanumesh + Mr. Varadareddy Name and Signature of Engineer : Hanumesh
 Name and Signature of User : D. PHANI KUMAR
27.11.19
 Mobile : 7093001023
 Email : hanumesh@srico-labworld.com Company Stamp

Scientific Research Instruments Company Private Limited

Telephone No : +91 9900055879 E-Mail : service@srico-labworld.com Website : www.srico-labworld.com

Branches: Bangalore: BBSR: Delhi: Hyderabad: Thane:
 Customer Copy : Head Office Copy: Branch Copy: Service Engineer Copy:
 Internal office use :