





## INSTALLATION REPORT

<b>Installation Date</b>		31-Oct-22		<b>Spincot ID</b>		S0001/12402		<b>Job Status (Completed/Pending)</b>		Completed	
<b>Instrument Model</b>		NA		<b>P.O.Number</b>		NA		<b>P.O Date</b>		NA	
<b>Principals</b>		NA		<b>Invoice No</b>		NA		<b>Invoice Date</b>		NA	
<b>Customer ID</b>				<b>Other Details :</b>		NA					
<b>Account Information</b>						<b>Contacts (Manager &amp; Users)</b>					
<b>Account Name:</b>		Indian Immunologicals Limited				<b>Designation</b>		<b>Mobile</b>		<b>E-mail</b>	
<b>Address &amp; City:</b>		Survery no-281-284&321, Biotech Park Phase 3, Karikapatla, Markook, Siddipet(DT), Telangana-502281				Mr.Kishan		Deputy Manager		8143260036	
<b>Phone :</b>											
<b>Mobile :</b>		8143260036									
<b>Email :</b>		k.kishan@indimmune.com									
<b>Website :</b>											
<b>Pre-Installation Checks</b>											
Line - Neutral (230 V ± 10 %, ± 1 V)		NA		Temperature (18 °C to 25 °C)						NA	
Line - Earth (230 V ± 10 %, ± 1 V)		NA		Humidity (40 % - 70 %)						NA	
Neutral - Earth (< 1 V)		NA		A.C (Yes/No)						Yes	
Instrument Chassis to Ground (< 100 Ω)		NA		Dust Free (Yes/No)						Yes	
<b>Other Observations, If any</b>											
<b>Installation &amp; Demonstration (Yes/No)</b>						<b>Qualification (Yes/No)</b>					
1. Installation Stacking		NO		1. Document Fill-up						Yes	
2. Flowline/Gasline Setup		NO		2. Document Verification						Yes	
3. Leak Pressure Check		NO		3. Is Verified Document Submitted						Yes	
4. Cables Organization & Tying		NO		4. Document submitted to						Mr.Kishan	
5. Solvent / Drain bottles - Capped		NO									
6. PC Settings		Yes									
7. PC Check		Yes									
8. Software Check		Yes									
9. Hardware Check		NO									
10. Demonstration		NO									
<b>Date &amp; Time</b>						<b>Instrument Configuration</b>					
<b>Date</b>		<b>Time In</b>		<b>Time Out</b>		<b>Model Name</b>		<b>Version</b>		<b>Serial Number</b>	
28-Oct-22		2:00 PM		5:30 PM		LABSOLUTIONS DB		6.11		1.54755903464	
31-Oct-22		10:30AM		6:00 PM							
<b>Engineer Observations :</b>						<b>Customer Comments :</b>					
Installation is Completed Properly						Software update done					
<b>Engineer Signature:</b>						<b>Customer Signature with Date &amp; Stamp :</b>					
											
<b>Engineer Name</b>		Mr.Saikumar.P				<b>Customer Name</b>		Mr.Kishan			
<b>Mobile</b>		7358773662				<b>Designation</b>		Deputy Manager			
<b>E-mail</b>		sotharaju.saikumar@spincotech.com				<b>Mobile</b>		8143260036		<b>E-Mail</b>	
										k.kishan@indimmune.com	

# License Certificate

## LabSolutions LC/GC

---

Upgrade from LCsolution

Serial Number : L54755903464

Security Code : JMY

This document certifies that you are legally permitted to use the software specified above ("SOFTWARE") subject to the terms and conditions of the attached SOFTWARE LICENSE AGREEMENT of Shimadzu Corporation ("SOFTWARE LICENSE AGREEMENT").

1 By installing or otherwise using SOFTWARE which accompanies this license, you agree to be bound by the terms and conditions of this document and SOFTWARE LICENSE AGREEMENT.

2 Total number of online systems shall not exceed the number set forth above. You agree to allow Shimadzu Corporation or its designee to audit the number of online systems upon reasonable notice during business hours.

Please keep this document in a safe place.

Shimadzu Corporation