



FirstSource Laboratory Solutions LLP
 Plot No.2, Block-D, Street No.3
 Kakatiyanagar, Hubsiguda, Hyderabad-500 007.
 Ph: 040-27175520-22 Fax: 040-27175523
 Email: sales@firstsourcels.com, Web: www.firstsourcels.com

INSTALLATION REPORT

CRM PM No:	Spinlab ID:	Engineer Name: <i>Jagan</i>	Date: <i>18/4/2022</i>
Account: <i>112</i>	Customer ID:	Letter	Date:
Address: <i>Gachibowli</i>	Principal: <i>Neelam</i>	Invoice no:	Date:
State: <i>TS</i> Pin Code:	Product Name: <i>Vortex</i>	Job Completed: Yes <input type="checkbox"/> No <input checked="" type="checkbox"/>	
SID In-Charge Name: <i>Narayana Pantu</i>	System User Details		
Designation: <i>Sr. Manager</i>	Name	Designation	Mobile
Department: <i>R&D</i>			
Phone:			
Fax:			
mobile:			
Email: <i>p.narayana@indimma.com</i>			
Website:			

Details: Installation *Neelam iSwix VT S/No. QJ049456*

Installation Procedure

- (i) Unpacked the instrument and checked physical conditions.
- (ii) Verified power connections and connected the power cord.
- (iii) Performed standard operation and instrument is ready to use.

AC Line Voltage Reading: P-N = 230 V P-E = 225 V E-N = 1 V UPS Stabiliser

Further Action Required Make permanent earth connection before start using the instrument.

This is to certify that the above mentioned work has been carried out and the system is working to our satisfaction.

Customer Comments:

Date			Time IN			Time OUT			Engineer Signature:		
<i>18/4/22</i>									<i>Jagan</i>		
									Customer Signature: (With Seal)		
									Date: <i>18-04-22</i>		