



TRIDENT EQUIPMENTS PVT. LTD.

A - 105,, Kailas Industrial Complex, Hiranandani Godrej Link Road, Vikhroli (W),
Mumbai - 79, Tel.: 2518 1704 / 05 / 06 Telefax : 2518 5463 E-mail : info@tridentequipments.com

Service Report

No. 63098

Name of Customer : <u>Tadras Instruments</u> <u>limited</u>	Product : <u>Thermo Scientific</u>
Address : <u>Cochin</u>	Model(s) : <u>omni star A215</u>
Contact Person/Dept : <u>Narajana</u>	Serial No.(s) : <u>X73452</u>
No. :	Installation <input checked="" type="checkbox"/> Warranty <input type="checkbox"/> Chargeable <input type="checkbox"/> Courtesy <input type="checkbox"/>
Phone : <u>9502033400</u>	Repaired at - Customer Site <input type="checkbox"/> Service Centre <input type="checkbox"/>
Email :	AMC : Preventive Maintenance Visit No. <u>-</u>
	Breakdown Visit No. <u>-</u>
	Date : <u>10.03.2023</u>

Nature of Problem :

Installation of omni star A215 pH/conductivity

Work Done :

- 1. Installed All Accessories as per instructions manual
- 2. Configured all parameters as per installation
- performed calibration with master calibration and with pH buffers
- And also calibration with resistance kit and with standards. observed results within limits
- Demo and training on instrument provided to user

Technician's Comment :

Installation of omni star A215 pH/conductivity successfully completed.

Customer's Comment :

Mr. N. Vinod Kumar
No. 01.02.2023
Engineer's Name & Signature :

U. Manikanta
10.02.2023
Customer's Name, Signature & Stamp :

8976898468