

SERVICE REPORT



G-315, Balaji Nagar, Kukatpally, Hyderabad - 500 072,
Ph (O) : 9505555984, E-mail : krishna@tecpharma.in

Date : 24/11/2020

S/R : 1611

Engineer's Name : D. Ashok

Mobile No. : 8897885376

Customer's Name :

Address Indian Immunological Limited
Wahapuram, Gachibowli post
Hyderabad - 500032
Telangana, India

Type of Services

- | | | | |
|-----------------|-------------------------------------|-----------------|-------------------------------------|
| 1. Installation | <input checked="" type="checkbox"/> | 6. Paid Service | <input type="checkbox"/> |
| 2. Breakdown | <input type="checkbox"/> | 7. Calibration | <input type="checkbox"/> |
| 3. AMC | <input type="checkbox"/> | 8. IQ/OQ/PQ | <input checked="" type="checkbox"/> |
| 4. Warranty | <input type="checkbox"/> | 9. Application | <input type="checkbox"/> |
| 5. Inspection | <input type="checkbox"/> | 10. Goodwill | <input type="checkbox"/> |

No.

Mail :

Job Completed : Yes No

Instrument

Model No.

Sr. No. / INST ID

Haier Ultra Low Temperature Freezer DW-86L729BP BEOGY8EBA00QGL880006/
B5/ULF/219

Problem As Reported by the Customer : Haier ULT Freezer installation

Service Done : Haier Ultra Low Temperature Freezer installation, completed, performed IQ, OQ and PQ qualification during 72 hours temperature mapping without load, 72 hours temperature mapping with load -70°C ± 1°C, power failure study and door open study done, now the Ultra Low Temperature Freezer working satisfactorily.

Spare Parts Used / Required

/
NA
/

To Be filled by Customer

- Instrument Installed & Found Working Satisfactorily
- Prevention Maintenance Carried Out Satisfactorily
- Fault Rectified & Instrument Working Satisfactorily
- Please await our written approval for parts replacement
- Calibration done Satisfactorily
- Comments (If Any)

Name : P. Subramaniam

Handwritten Signature

Date : 24-11-20

Customer Signature

Engineer's Signature : D. Ashok
24/11/2020

* Only for parts Changed & all other charges will be extra as applicable

Parts replaced are chargeable except during warranty period. Consumables used for the instrument are not covered by warranty & hence are chargeable. Parts replaced due to mishandling during operation or otherwise will also be charged in case. In case of any difficulty or suggestion, please contact Service Manager.