

Installation report

Client: M/s Indian Immunologicals Limited

Date: 28-Feb-2022

Indian Immunologicals Limited called DrM technical persons for installation of FUNDAMIX® vibromixer type serial number 3172

Frist 1000L water was filled in the fermenter (E.Q ID: B8/FER/085) & Vibromixer was installed. Shaft distance is adjusted 4 mm between the motor & shaft. FUNDAMIX® Size 3 ran to check the efficiency as per the test report FO 3.3.8 maximum amplitude 2.4 mm is achieved at 230 volts.

Performance of FUNDAMIX® FM-3 checked after installation.

Potentiometer	Voltage (V)	Amplitude (mm)
1	137.8	1.0
2	148.8	1.15
3	161.7	1.30
4	174.9	1.45
5	192.0	1.80
6	208.5	2.00
7	222.7	2.15
8	230.0	2.30
9	235.6	2.40

EVR Reddy
28.02.22


28.02.22



